

Maestro 3D Ortho Studio Viewer Specification Requirements

The Maestro 3D Ortho Studio Viewer software is designed to operate with the following computer hardware/software requirements:

PC : - A CPU with 2.0Ghz or greater (an Intel Core 2 Duo)
- Minimum of 2 GB RAM

Operating System: Windows XP, Windows Vista or Windows 7.

Video Card: OpenGL compatible. For best results, we recommend using a Nvidia video card. AVI and integrated video cards are not recommended as they do not support all OpenGL specifications.